

## **MicroFit HealthStar Manager - Version 5.7.1**

**May 2011**

### **HealthStar Manager 5**

#### Export to District Database Enabled

In Administration > Database there is now an Export button. This feature will send (export) student fitness data from the school database up to the district database located on MicroFit servers. After all schools have completed the export, the district administrator may use HealthStar Manager 6 software to view and print summarized fitness scores for all students in the district.

Note: HealthStar Manager 6 is a different software package from HealthStar Manager 5. It is designed for district fitness reporting only. Contact MicroFit customer service (or technical support) for information on HealthStar Manager 6.